

RespiClick® Inhaler

Medicines: ProAir® (albuterol sulfate), Air Duo® (fluticasone/salmeterol)



ProAir® (albuterol sulfate)



Air Duo® (fluticasone/salmeterol)

The RespiClick is an inhaler that puts a dry powdered medicine into your lungs.

- **Do not use with a spacer or valved holding chamber.**
- **Do not shake** the RespiClick
- **Do not prime** the RespiClick before using it.
- Make sure the cap is closed before and after each dose.
- Don't open the cap until you are ready to take your dose.

How to use the RespiClick inhaler

1. Check the counter on the inhaler before use
2. Open the cap.



3. Inhale



- Hold the inhaler upright.
- **Open the cap all the way until you hear a “CLICK.”**
- Your RespiClick is now ready to use.
- Face away from the RespiClick and breath out fully. **Do not breath out into the inhaler mouthpiece.**
- **Don't let your lips or fingers block the vent above the mouthpiece.**
- Take a quick deep and full breath in.
- Hold your breath for 10 seconds.
- Remove the inhaler from your mouth and check the counter on the back of the inhaler to make sure you got the dose.

4. Close the cap.




- **Always close the cap** after each dose so your inhaler will be ready for your next one.
- Repeat these steps if more than 1 puff is recommended by your doctor.

How to clean your RespiClick

- **Do not wash or put any part of the RespiClick in water.**
- Clean mouthpiece weekly with a clean, dry tissue and cotton swab.

More Information

- If the cap breaks off, the inhaler can no longer be used and you will need a new one.
- Store it in a dry area, at room temperature, with the cap closed. Don't keep it anywhere that is too hot or too cold.
- When using a dry powder inhaler, you may not taste the powder and you may not feel the medicine being delivered.
- There is a counter on the back of the inhaler. When the counter says 20, you should call the pharmacy for a refill.
- When the counter says 0, the RespiClick is empty and should be thrown away and replaced with a new RespiClick



To watch a video on how to use the RespiClick® Inhaler, scan the QR code below.



<https://bit.ly/3UqwVAF>

childrenscolorado.org

Author: Asthma Education | Approved by Patient Education Committee | Valid through 2026

The information presented is intended for educational purposes only. It is not intended to take the place of your personal doctor's advice and is not intended to diagnose, treat, cure or prevent any disease. The information should not be used in place of a visit, call or consultation or advice of your doctor or other health care provider.