

Name:	DOB:	School:	Grade
Parent/Guardian:	Phone Number/s:	Provider:	Phone:

Pump Bolus Feed:

Supplies:

1. Gloves
2. Formula
3. Pump and feeding bag
4. G-tube extension set
5. ____ mL syringe for water flush, if ordered.

Procedure:

1. Place Supplies on clean surface.
2. Wash hands and put on gloves.
3. Observe for distention. If the stomach looks larger than normal, call the parent and RN for further instructions.
4. Pour ____ mL formula into feeding bag.
5. Attach g-tube extension set to feed tubing.
6. Insert tubing into pump as instructed by RN and turn pump on.
7. Prime tubing with formula, using the prime button on the pump, stopping at the end of the g-tube extension set. Ensure there are no bubbles in the line.
8. Ensure pump is set at ____ rate.
9. Unclamp the g-tube extension set and feeding tube and push start on pump.
10. Once feed is complete disconnect feed tubing from extension set and flush g-tube extension set with ____ mL of water, if ordered.
11. Disconnect g-tube extension from g-tube and close safety plug.
12. Remove gloves and wash hands.
13. Rinse feeding bag and extension tube thoroughly with warm water and allow to air dry. Warm, soapy water may be used if necessary (e.g. visibly soiled) but should be rinsed thoroughly after. Store formula as instructed.
14. Document feed on G-tube feeding log.

Additional Comments/Instructions:

This document, and the information it contains was created by Children's Hospital Colorado ("CHCO") to serve as a guideline and reference tool for use by CHCO employees while acting within the scope of their employment with CHCO. The information presented is intended for informational and educational purposes only. It is not intended to take the place of your personal physician's advice and is not intended to diagnose, treat, cure or prevent any disease. The information should not be used in place of a visit, call, consultation or advice of your physician or other health care providers.

