

Name:	DOB:	School:	Grade
Parent/Guardian:	Phone Number/s:	Provider:	Phone:

Syringe Push Bolus Feed:

Supplies:

1. Gloves
2. Formula
3. 60mL syringe
4. G-tube extension set
5. ___mL syringe for water flush, if ordered.

Procedure:

1. Place supplies on clean surface
2. Wash hands and put on gloves
3. Observe for distention. If the stomach looks larger than normal, call the parent and RN for further instructions.
4. Ensure the g-tube extension set is flushed with water.
5. Attach extension set to the g-tube
6. Draw up formula into the 60mL syringe and connect it to the feeding extension set.
7. Unclamp the g-tube extension set and push ___ mLs over ___ minutes.
8. Clamp extension set and draw up the rest of the formula into the 60mL syringe and repeat steps 7 and 8 until feed is complete.
9. Flush g-tube extension with ___ mL of water after feeding is complete, if ordered.
10. Disconnect g-tube extension from g-tube and close safety plug.
11. Remove gloves and wash hands.
12. Rinse 60mL syringe and extension tube thoroughly with warm water and allow to air dry. Warm, soapy water may be used if necessary (e.g. visibly soiled) but should be rinsed thoroughly after. Store formula as instructed.
13. Document feed on g-tube feeding log.

Additional Comments/Instructions: